


iPro Quick Reference Guide


1. Enter username and password. Turn on **Remember Me** to save your username




2. Tap the Today's Cases tab to view all the cases.
 - A. Search for a case in the search field
 - B. Change the facility to SUNY Downstate by selecting the dropdown option at the top right



3. Select your case
4. Then tap **Take Case**. This will assign you to the patient's chart as a "Participant"



5. Begin charting
6. Tap **Checklist** from the left side menu and start charting



Isabel Phillips | 257142 | 80 Yrs F | 12/31/1937 | 75 kg

Revision of femoral component | Penicillin

Physical Status: ASA 2 | Anes Start: [] | Time-Out: [] | Incision: [] | Surgery End: []

Pre-Procedure

PT Double Identified H&P Reviewed Chart Reviewed Consent Signed

Patient Re-Evaluated Prior to Induction Pre-Induction Vitals are First Set of Vitals

Time Out Performed Prior to Induction

Pre-Anesthesia State

Awake Asleep Apprehensive Unresponsive

Patient Safety


Eye Care

Ointment Taped Goggles Closed w/o pressure

Position

-- Select -- Pressure Points Checked/Padded Axillary Roll Shoulder Roll

7. Add **Blocks** if needed in the **Regional** section
 - A. User can add additional blocks by selecting the + symbol in the right-hand corner of the Block section
 - B. User can take pictures of Ultrasound Guidance with the camera icon and save it with the patient record



PROVIDERS

PREOP

CHECKLIST

REGIONAL

CHART

EMERGENCY

MACRA

NOTES

SIGNATURES

POSTOP

AQI 53 - Documentation of Anticoagulant and Antiplatelet Medications when Performing Neuraxial Anesthesia/Analgesia or Pain Procedures

Patient on anticoagulant/antiplatelet Tx Not on anticoagulant/antiplatelet Tx

7 Block/Neuraxial Block **A** +

Block	Site	Approach	Position
Femoral Block	Left	Anterior	Supine

Type	Gauge	Catheter at skin(cm)	Attempts
Tuohy	16	3	1

Epidural needle at [] cm by LOR technique Air Saline


Local Anesthetic []

Test Dose [] cc Results [] at dose Negative Positive


Ultrasound Guidance **B** [Camera Icon] Sterile gel and sterile probe covers used

Procedure Note

8. To begin charting Vitals, Gases, Monitors, Fluids, and Vent Settings open the **Chart** section
 - A. Enter your **Anes Start Time**
 - B. To add Gases, Drugs, Monitors, Vent Settings, or Fluids tap the **+** symbol on the grey header of the item you wish to add and select from the dropdown of available options



9. To start an **Auto Vitals**, tap the Auto-Vital button
 - A. Scan the barcode on your anesthesia machine
 - B. Or tap **List** to get a list of machines
 - C. Turn on or off any additional readings that you wish to document



10. Tap **Finalize** to find if the user is missing any required items

