



CERNER BRIDGE SPECIMEN COLLECTION

BAR-CODE SCANNING UTILIZING HAND-HELD ZEBRA DEVICE

INSTRUCTIONS

1. Assemble the following before entering the patient's room:
 - a. Specimen Container (blood tubes, specimen container)
 - b. Zebra Hand-held Device
 - c. Printer



2. Use the devices to scan the patient ID and **ASK** patient to verify their **FULL NAME** and **DATE OF BIRTH**.
3. Proceed with specimen collection (steps outlined in the video below)

LINK TO VIDEO:

[Watch my Powtoon: Video: Collecting Specimen Cerner Bridge](#)

QR CODE TO VIDEO:



ATTESTATION: Attest you have viewed this educational video by completing this short form:

<https://forms.gle/fDqPkaiJKJMAhSdd7>



ATTESTATION QR CODE