

Gastroenterology and Advanced Endoscopy Phone (718) 270-4772 • Fax (718) 270-7201 www.DownstateGl.org

Patient Instructions for PillCam Small Bowel Capsule Endoscopy with SensorBelt

	☐ Monday	☐ Tuesday	☐ Wednesday	☐ Thursday	☐ Friday
Date:			_	Time:	
			LOCATION:		
		Univers 470 Cla Brookly (718) 2 Suite A	Downstate Central sity Hospital of Bro arkson Avenue yn, NY 11203 170-4788 Copy Center - Suite	oklyn	
		339 Hid Brookly (718) 7 Admitt Subway	Downstate Long Is cks Street (corner of yn, NY 11201 780-1991 ing Office (then pro y: 2,3,4,5,J,M,R to I B61 or B63	of Atlantic Ave.)	or—Endoscopy Suite)
		9036 7 Brookly (718) 5 Urgent Subway Buses:	Downstate Bay Rid th Avenue (at 92nd yn, NY 11228 667-1464 Care Center (then y: R to 95th Street B70 to 92nd Street B16 to 95th Street	d Street) proceed to 6th F	

Patient Instructions for PillCam Small Bowel Capsule Endoscopy with SensorBelt

A PillCam small bowel capsule endoscopy enables your doctor to examine your small intestine. You will swallow a vitamin-sized capsule that contains a miniature camera. While the video capsule travels through your body, images are sent to a data recorder that you wear on a belt around your waist. The capsule endoscopy will last approximately 8-9 hours. Afterwards, your doctor will view the images on a video monitor.

1 Week Before Your Exam:

• Discontinue medications containing iron. This includes multivitamins with iron.

Day Before Your Exam:

At noon, begin a clear liquid diet (no red, purple or dark-colored liquids):

- Water
- Clear broth or bouillon
- Coffee or tea (without milk or non-dairy creamer)
- Gatorade, Pedialyte or Powerade
- · Carbonated & non-carbonated soft drinks
- · Kool-Aid or other fruit-flavored drinks
- Strained fruit juices (no pulp)
- JELL-O, popsicles, Italian ices, ice pops and hard candy

10 Hours Before Your Exam:

• Stop drinking all liquids, except for necessary medication with a sip of water.

Day of Your Exam:

- You may take medications up to 2 hours before your exam with a sip of water.
- We will be placing a belt with sensors across your abdomen. Please wear an upper garment of thin, natural fiber (such as a T-shirt) that is long enough to reach at least hip level.

After Swallowing the PillCam Small Bowel Capsule:

- Be sure the SensorBelt is tight at your waist. Do not attach anything to it.
- You may drink colorless liquids starting 2 hours after swallowing the capsule.
- You may eat a light snack 4 hours after capsule ingestion. After the exam is completed you may return to your normal diet.
- Every 15 minutes Check the blue flashing DataRecorder light to be sure it is blinking twice per second. If it stops blinking or changes color, note the time and contact your doctor.
- Do not disconnect the equipment or completely remove the DataRecorder at any time during the procedure.
- Treat the DataRecorder with care. Avoid sudden movements or banging of the DataRecorder.
- After ingesting the capsule and until it is excreted, you should not have a Magnetic Resonance Imaging (MRI) examination or be near an MRI device.
- The capsule is disposable and passes naturally with your bowel movement. You should not feel any pain or
 discomfort during capsule excretion. (If you are not sure that the capsule has passed out of your body and you
 develop unexplained nausea, abdominal pain or vomiting, contact your doctor.)

After Completion of the Exam:

- Return to the doctor's office to have the equipment removed.
- Or, if instructed to remove the equipment yourself:
 - Remove the SensorBelt and DataRecorder and keep in a safe place.
 - Handle the equipment carefully. Do not expose it to shock, vibration or direct sunlight.
 - Return all equipment to your doctor's office as soon as possible.
- The images acquired during your exam will be downloaded to a computer for physician review. After the doctor has looked at the video, you will be contacted with the results.

Time to return to facility:	Special Instructions:
Emergency contact:	